The background of the slide is a high-quality photograph of various surgical instruments, including several pairs of forceps and a scalpel, arranged on a light-colored surface. The instruments are made of polished metal and are slightly out of focus, creating a sense of depth. A dark red circular graphic is overlaid on the left side of the image, partially obscuring some of the instruments.

Anterolateral Approach

Instrumentation rationale

Introduction

To maximize functional recovery, the minimally invasive, Anterolateral Watson-Jones approach was developed to minimize the postoperative dislocation rate while maintaining the abductor strength.

It uses the classic Watson-Jones interval, but without detaching, cutting, or otherwise damaging the abductor. In order to help surgeons better utilize this approach, MicroPort Orthopedics has developed

specialized instrumentation to take advantage of the surgical technique and maximize the operative window. When used in conjunction with intraoperative fluoroscopy and MicroPort's product portfolio, these instruments optimize the positioning of MicroPort Orthopedics' implants to deliver better results.



Optimized surgical window

The MicroPort Orthopedics Anterolateral instruments are designed to maximize a surgeon's visibility in a smaller surgical window. Approach-specific retractors increase the visible operative field and protect soft tissue to aid in recovery. Low profile, offset handles facilitate the proper preparation of the patient's bones and implantation of the final implants without putting extra strain on the patient's position.

Comprehensive instrumentation

This instrumentation system can be adapted to match an individual surgeon's needs and technique. The variety of retractors, offset handles, and preparation instruments provide a comprehensive system that can be tailored to the experience and preferences of the surgeon.

Minimally invasive implant portfolio

As a leader in minimally-invasive surgical techniques MicroPort Orthopedics has designed these instruments to work in tandem with our established product portfolio. The Profemur® stem line offers implants such as the Profemur® Preserve and Profemur® Gladiator that have been designed to provide great results in a smaller operative window. The Prime® acetabular system offers intraoperative flexibility to deliver optimal results in a wide range of patients.

Full Function *Faster*®





Optimized surgical window

Increased exposure

Specialized retractors are designed to utilize anatomic structures to protect soft tissue and increase the surgical window.

Right Angle Cobra Retractor, Curved Cobra Retractor, Double Prong Retractor, Short Hohmann Retractor, and Arch Retractor.



Reproducible preparation

Simple instrumentation

Offset broach handles take advantage of the patients positioning and allow the surgeon to prepare the femur without placing additional strain on the soft tissues.


- 1.4" lateral offset from broach to handle
- 1" and 2" anterior offset options
- 20° angle offset from broach handle to clear body and correct potential varus
- Strike plate centered on handle to minimize eccentric force on femur
- Extraction plate separate to facilitate broach removal
- Modular attachment point for inserter, broaches, and box chisel



Comprehensive instrumentation

Efficient tray design




Supplemented by the desired implant-specific instrumentation, only three Anterolateral instrument trays are needed for each case.



Core instruments

Retractors






Broach handles



Broach handle options

MicroPort offers straight, 1" offset angled, 2" offset angled, and double offset broach handles to accommodate surgeon preference and optimize the surgical technique

- 1" offset angled broach handle
- 2" offset angled broach handle
- Double offset angled broach handle




Broach handle comparison

Features	1" offset angled	2" offset angled	Double offset
Anterior offset	20°	20°	1: 45° 2: 2°
Lateral offset	2.0"	2.0"	2.0"
Overall length	10.3"	10.3"	12.0"
Extraction strike	Tap-out hump	-	Underside of plate


Offset cup impactor

The offset cup impactor with a modular version guide allows the surgeon to adjust the insertion angle for each patient, while maintaining the appropriate version during implantation for reproducible results.



Offset cup impactor with quick release


The offset cup impactor with a quick release attachment allows for quick release of the handle with one simple quarter turn following cup impaction. This cup impactor also allows for the attachment of a modular version guide.




Minimally invasive product portfolio

Profemur® Preserve

Featuring 24 stem sizes, optimized neck geometries to maximize head center coverage, and a reduced lateral shoulder, this short stem is well-suited for minimally invasive procedures. The curved medial geometry facilitates load transfer to the cortical bone and helps to ease insertion, and direct lateralization of head centers allows intraoperative offset adjustments to be isolated.



Straight
34mm
Varus 8°
37mm



Straight
31mm
Varus 8°
33mm



Straight
27mm
Varus 8°
29mm





■ **Size 1-4**
■ **Size 5-8**
■ **Size 9-12**

The difference in neck length between each group is roughly one head size, 3.5 mm.

The clinical centroid is the driving force behind the Preserve design. A study of 900 radiographs, sampled globally, confirmed that femur size and neck length are positively correlated, and revealed that three neck groups is the optimal number for minimizing implant options while maximizing head center coverage and allowing for improved operating room flow. The design was robustly validated with an additional 150-case randomized templating study.

Profemur® Gladiator

The triple-tapered wedge shape and the reduced lateral shoulder on the Profemur® Gladiator make it a great option for minimally-invasive procedures such as Direct Anterior. Additionally, horizontal grooves are designed to more evenly distribute load forces, while vertical grooves assist with rotational stability.



Profemur® TL

Designed with a dual-taper philosophy, the Profemur® TL stem features a slim anteroposterior dimension that preserves bone and provides stable, three point fixation in the proximal femur. The reduced lateral shoulder, narrow profile, and broach only technique combine to make this stem a great option for minimally invasive techniques such as Direct Anterior Total Hip Arthroplasty.

MicroPort offers implants to
match your **stem philosophy**

View our full product line at microportortho.com



Liner options: standard, lipped, and
lateralized/face changing

A-Class: highly cross-linked
E-Class: vitamin E blended polyethylene

Prime Acetabular System


The Prime Acetabular Cup System is the next step in the evolution of the successful Dynasty® Acetabular Cup System. It is the first system that supports a variety of surgical approaches and is optimized for a highly cross-linked polyethylene bearing surface, eliminating the compromises associated with alternative bearings.

BioFoam® Cancellous Titanium coating

Product ordering information

ATALKIT1 ANTEROLATERAL RETRACTOR KIT

- A. Right Angle Cobra Retractor (3△)
20162000
- B. Short Hohmann Retractor (5△)
20162001
- C. Curved Cobra Retractor (6△)
20162004
- D. Arch Retractor (0△)
20162006
- E. Double Prong Retractor (4△)
20162010



ATRRKIT2 ANTERIOR 2" OFFSET BROACH HANDLE KIT

- A. 2" Offset Angled Broach Handle–Left
ATBRHANL
- B. 2" Offset Angled Broach Handle–Right
ATBRHANR



A collection of orthopedic surgical instruments, including a Muller neck elevator, a bone hook, a starter rasp, and a T-handle, arranged on a perforated metal tray. The instruments are labeled with letters A through K.

-

M. Modular offset stem impactor
PRMOD460



Full Function, Faster®



MicroPort Orthopedics Inc.
5677 Airline Road
Arlington, TN USA 38002
866 872 0211

microportortho.com

The CE-Marking of Conformity is applied
per catalog number and appears on the
outer package label, if applicable.

Trademarks and Registered marks of
MicroPort Orthopedics Inc.
© 2019 MicroPort Orthopedics Inc.
All Rights Reserved. 016408B AUG2019