



HIP PRODUCT BROCHURE

Table of Contents

CUPS	1	Prime™ Acetabular Cup System
	2	Dynasty® Acetabular Cup System
	2	Head Options
	3	Gladiator® Bipolar Acetabular System
STEMS	4	Profemur® Preserve Stems
	5	Profemur® TL2 Stem
	6	Profemur® TL Stem
	7	Profemur® Gladiator® Stems
	8	Profemur® Z and Plasma Z Stems
	9	Profemur® Renaissance® Classic Stem
	10	Profemur® R Stems
MATERIALS	11	BioFoam® Cancellous Titanium
	13	A-Class® Highly Cross-linked Polyethylene
	14	E-Class™ Vitamin E Blended Highly Cross-linked Polyethylene

Prime™ Acetabular Cup System

The Prime™ acetabular cup system is the next step in the evolution of the successful Dynasty® acetabular cup system. The system is optimized for a highly cross-linked ultra-high molecular weight polyethylene bearing surface, eliminating the compromises associated with modularity to accommodate alternative bearing surfaces. By focusing on a singular bearing surface, the shell has also been optimized for modern bone fixation surfaces. Designed with simple, versatile instrumentation, the system can be used with a variety of surgical approaches.

DESIGN FEATURES

PRIME™ SYSTEM	
Shell sizes	46 – 68mm in 2mm increments
Liner diameter	22, 28, 32, 36, 40, and 44mm
Bearing surfaces	A-Class® Highly Cross-linked Polyethylene and E-Class™ Vitamin E Blended Highly Cross-linked Polyethylene
Shell options	Solid, 3-hole MH
Shell coating	BioFoam® Cancellous Titanium
Liner types	Standard, lipped, and lateralized/face changing



Robust Locking Mechanism

Allows for multiple liner configuration options, promotes easy liner insertion, and maximizes push-out strength

Optimized Screw Hole Locations

Create divergent fixation and allow for 18° of screw angulation

Square Impaction Dimple

Provides rotational control during implantation and functions with a quick release impactor mechanism

Optimized Head to Shell Ratio

Allows the use of a 36mm head and liner in a 50mm shell with no compromise to liner thickness

Minimized Shell Thickness

Decreases stiffness, discourages stress-shielding, and allows for optimized liner thickness

Dynasty® Acetabular Cup System

The Dynasty® acetabular cup system offers the benefits of cementless press-fit design together with its clinical established A-Class® highly cross-linked polyethylene. The configuration of Dynasty® shells allows for up to 10 screw holes on the larger diameter cups in conjunction with the latest bone fixation technology, BioFoam® cancellous titanium, making this cup the ultimate primary and revision option.

DESIGN FEATURES

	DYNASTY® POROUS-COATED SYSTEM	DYNASTY® BIOFOAM® SYSTEM
Shell sizes	46-68 in 2mm increments	46-76 in 2mm increments
Poly liner diameters	28-46	28-54
Bearing surface	A-Class® Highly Cross-linked Polyethylene	A-Class® Highly Cross-linked Polyethylene
Number of screw holes	3	3, 7, 8, 10 depending on cup size
Shell coating	Porous beads	BioFoam® Cancellous Titanium
Revision poly option	36mm ID from a 52-68mm	36mm ID from a 52-68mm
Liner types	Standard, 15 degree	Standard, 15 degree



HEAD OPTIONS



CoCr metal



BioloX® Delta®
ceramic head



BioloX® Delta®
option head



A-Class®
backfilled metal



BCH®
ceramic head



Titanium Sleeve
used with BioloX® Delta®
option head



Titanium Sleeve
used with Backfilled metal
& BCH® ceramic heads


Gladiator® Bipolar Acetabular System

Gladiator® bipolar acetabular system is a bipolar hip implant design that features a cross-linked polyethylene bearing surface with an enhanced lock detail designed for strength. Historical concerns with traditional bipolar designs have included loosening of the insert, disassociation of the head from the shell, and osteolysis resulting from polyethylene wear. This system is designed to address these concerns to give surgeons greater confidence when using a bipolar implant.

There is an UHMWPE support ring inside the shell that is permanently fixed. There is also an UHMWPE locking ring that assembles above the support ring and locks into place once the head is inserted into the shell.

DESIGN FEATURES

GLADIATOR® BIPOLAR ACETABULAR SYSTEM	
Shell sizes	36 – 59mm
Head diameters	22 – 36mm
Bearing surface	A-Class® Highly Cross-linked Polyethylene



Design features of Gladiator® Bipolar Acetabular System



Locking ring in its "unlocked" position, without the head in the shell.



Head has been inserted. This causes the locking ring to engage and locks the head into place.

Profemur® Preserve Stem

DESIGN FEATURES

ABBREVIATED TECHNIQUE: BROACH ONLY

Broach to templated size

Implant size corresponding to broach



SIZES

Neck	Optimized length
Stem	1-12

For additional risk information, please consult the Instructions For Use package insert.

DESIGNED TO MAXIMIZE HEAD CENTER COVERAGE

The design team reviewed over 900 radiographs to determine Preserve modular stem size, modular neck type, and head size in an effort to optimize component placement with existing patient anatomy. The templated neck/head combinations were converted into a frequency distribution to calculate the clinical centroid for each stem size.



Profemur[®] TL2 Stem

DESIGN FEATURES

ABBREVIATED TECHNIQUE: BROACH ONLY

Broach to templated size


Implant size corresponding to broach

SIZES

Neck	Optimized length for each size grouping
Stem	0–12

Available in Classic only

For additional risk information, please consult the Instructions For Use package insert.



Profemur[®] TL Stem

DESIGN FEATURES

ABBREVIATED TECHNIQUE: BROACH ONLY

Broach to templated size


Implant size corresponding to broach

SIZES

Neck Short and long lengths

Stem 1–12

For additional risk information, please consult the Instructions For Use package insert.



Profemur® Gladiator® Stems

DESIGN FEATURES

ABBREVIATED TECHNIQUE: BROACH ONLY

Broach to templated size

Implant size corresponding to broach

SIZES

Neck Medium length



Stem 1-10

PRESS FIT STEM SIZES

Neck Medium length

Stem 4,6,8,10,12

For additional risk information, please consult the Instructions For Use package insert.



Profemur[®] Z and Plasma Z Stems

DESIGN FEATURES

ABBREVIATED TECHNIQUE: BROACH ONLY


Broach to templated size

Implant size corresponding to broach

SIZES

Neck	Short and long lengths
Stem	1-9

For additional risk information, please consult the Instructions For Use package insert.



Profemur® Renaissance® Stem

DESIGN FEATURES

ABBREVIATED TECHNIQUE: REAM & BROACH

Ream to templated size or cortical chatter

Sequentially broach with reduced flare broaches to corresponding reamer size

If it is not secure, switch to standard flare broach or ream and broach to larger stem size

Implant size and flare corresponding to broach size and flare

SIZES

Neck	Short and long lengths
Stem	10-16 (Reduced flare) 10-18 (Standard flare)

For additional risk information, please consult the Instructions For Use package insert.

Neck options


Short and long neck lengths, straight (135° CCD) and varus 8° (127° CCD) neck angles allowing for multiple head center positions to meet range of anatomical needs

Flare

Standard and reduced flare options designed to maximize metaphyseal fit


Proximal taper

3° proximal anterior/posterior taper designed for proximal fill



Distal slot

Minimizes stem stiffness to discourage fracture during stem insertion



Driving platform

Dimple designed for unidirectional loading during stem insertion and oval slot designed for rotational control during stem insertion

Lateral shoulder

Rounded lateral shoulder designed to ease stem insertion and minimize risk of fracture

Plasma spray

Coating thickness provides 1.0mm (0.5mm per side) additional press-fit

Surface roughness

Titanium stem surface has glass-beaded texture

Distal splines

Designed to provide additional 1mm press-fit (0.5mm per side) for rotational stability

Distal bullet tip

Round distal tip designed to reduce the risk of fracture during insertion and minimize point contact after implantation

BioFoam® Cancellous Titanium


STRUCTURE

Cancellous or trabecular bone is composed of a network of rod- and plate-like elements that provide porous macrostructure for blood vessels and marrow.


The structure of Biofoam® cancellous titanium metal resembles that of trabecular bone. Made from commercially pure titanium, the pore cell size averages 530µm and the diameter of interconnecting pores averages 200µm. The porosity is between 60 and 70%, creating an open cell structure that allows deep bone fixation for long-term stability.




**BioFoam® Cancellous Titanium on Prime™
Acetabular Cup System**



Structure of BioFoam®



**Trabecular structure of BioFoam® Cancellous
Titanium, 200x original magnification**



**Trabecular structure of BioFoam® Cancellous
Titanium, 100x original magnification**


BioFoam[®] Cancellous Titanium

COMPRESSIVE STRENGTH & MODULUS


Compressive strength measures the maximum amount of compressive load a material can bear prior to fracturing.

Compressive modulus is a normalized measure of a material's stiffness measuring how much a material compresses under load without permanently deforming.

BioFoam[®] cancellous titanium is engineered to have a modulus similar to that of bone facilitating even and consistent bone loading to prevent stress shielding, promote long-term interdigitation, and enhance stability.



Compressive strength of Biofoam[®] metal as compared to bone and Trabecular Metal™




Compressive modulus of Biofoam[®] metal as compared to bone and Trabecular Metal™

FRICTIONAL PROPERTIES

Frictional resistance provides initial stability between the implant and bone, thereby helping prevent the implant from movement immediately following implantation. Immediate rigid fixation is crucial to the ingrowth process. If the implant is moving, it is no longer working as a stable construct for bone growth, and the bone will not be able to attach itself to the implant for long-term fixation.

In a study, BioFoam[®] cancellous titanium had a significantly higher coefficient of friction than porous tantalum (Trabecular Metal), plasma spray, and sintered beads (p=0.007, 0.051 and 0.001, respectively.)




Frictional properties of Biofoam[®] metal compared to traditional ingrowth coatings and Trabecular Metal™

A-Class® Highly Cross-linked Polyethylene


DESIGN FEATURES

- **No oxidation¹**
- **Undetectable free radicals²**
- **92% or greater reduction in wear³**



Wear rate per million cycles of A-Class® highly Cross-linked Polyethylene

MicroPort Orthopedics utilizes the following manufacturing processes for cross-linking its A-Class® polyethylene.



Poly material selection

Process begins with compression molded GUR 1020, which has a higher impact strength, tensile strength, and yield strength than GUR 1050⁷.

Heat treatment

Following irradiation, the rods are heated above the melting point of the polyethylene to eliminate residual free radicals, form additional cross-links, and improve the oxidative stability of the material⁸.

Cross-linking process

GUR 1020 rods are gamma irradiated to a dose of 7.5 MRads to facilitate cross-linking and enhanced wear resistance, but also maintain mechanical properties of the material.

Machining and final sterilization

Liners are machined, cleaned, packaged, and sterilized using ETO sterilization, which does not reintroduce free radicals or cause any other measurable change to the polymer.

E-Class™ Vitamin E Blended Highly Cross-linked Polyethylene

Joint replacement is occurring in increasingly younger patients. Over 20% of primary total hip arthroplasty (THA) procedures in the United States occur in patients 55 years of age or younger. This younger patient population expects to remain active throughout their lives, creating a need for longer-lasting and higher performance implants.


Reference: 2018 National Inpatient Sample, Hospital Cost and Utilization Project, Agency for Healthcare Research and Quality, US DHHS

DESIGN FEATURES

- Improves wear resistance⁵
- Undetectable free radicals⁶
- Maintains mechanical strength⁶



MicroPort Orthopedics utilizes the following manufacturing processes for cross-linking its E-Class™ Vitamin E blended highly cross-linked polyethylene.



Poly material selection

Process begins with compression molded GUR 1020, which has a higher impact strength, tensile strength, and yield strength than GUR 1050⁷.

Vitamin E treatment

Following irradiation, Vitamin E stabilizes any remaining free radicals and continuously prevents oxidation. The presence of Vitamin E eliminates the need for remelting, resulting in improved mechanical strength.

Cross-linking process

GUR 1020 rods are gamma irradiated to a dose of 10 MRads to facilitate cross-linking and enhanced wear resistance, but also maintain mechanical properties of the material.

Machining and final sterilization

Liners are machined, cleaned, packaged, and sterilized using ETO sterilization, which does not reintroduce free radicals or cause any other measurable change to the polymer.

1. Benchtop data on file at MicroPort Orthopedics.
2. Benchtop data on file at MicroPort Orthopedics.
3. Compared to conventional poly, Benchtop data on file at MicroPort Orthopedics.
4. Reference: 2018 National Inpatient Sample, Hospital Cost and Utilization Project, Agency for Healthcare Research and Quality, US DHHS
5. Compared to conventional polyethylene. Benchtop data on file at MicroPort Orthopedics.
6. Benchtop data on file at MicroPort Orthopedics.
7. "Standard Specification for Ultra-High-Molecular Weight Polyethylene Powder and Fabricated Form for Surgical Implants". ASTM International, 2014. F648-14.
8. Benchtop data on file at MicroPort Orthopedics.



MicroPort Orthopedics Inc.
5677 Airline Road
Arlington, TN USA 38002
866 872 0211

microportortho.com

The CE-Marking of Conformity is applied per catalog number and appears on the outer package label, if applicable.

Trademarks and Registered marks of MicroPort Orthopedics Inc.
© 2023 MicroPort Orthopedics Inc. All Rights Reserved.
011082E DEC2023