



# PRIME<sup>®</sup>

## ACETABULAR CUP SYSTEM




DESIGN RATIONALE

# Introduction

The Prime Acetabular Cup System is the next step in the evolution of the successful Dynasty® Acetabular Cup System. It is the first system that supports a variety of surgical approaches and is optimized for a highly cross-linked polyethylene bearing surface, eliminating the compromises associated with alternative bearings.

**The Prime Acetabular Cup System is the fifth generation of MicroPort Orthopedics' acetabular cup systems and is built on over 25 years of clinical heritage.**

**Magnification of trabecular structure of BioFoam® Cancellous Titanium**



BioFoam® cancellous titanium newest generation of cementless fixation

- The structure of the BioFoam® Cancellous Titanium acts as a biological scaffold designed to enhance bone apposition to support long-term implant strength and stability.<sup>1</sup>
- Porosity between 60–70% mimics natural trabecular architecture.
- Pore cell size averages 530µm and diameter of interconnecting pores averages 200µm providing a macrostructure for vascularization and bone apposition.
- Coarse titanium struts bite into bone providing initial scratch-fit and enhancing immediate rigid fixation.
- Compressive modulus will withstand impaction and repeated load and potentially transfer applied load to native bone, discouraging stress-shielding.<sup>2,3</sup>
- High abrasion resistance to withstand surgical impaction while still maintaining a 70% porous structure to allow osseous interdigitation.<sup>4</sup>



**Radiused shell edges**

Minimize soft-tissue irritation

**Robust locking mechanism**

Allows for multiple liner configuration options, promotes easy liner insertion, and maximizes push-out strength

**Optimized screw hole locations**

Create divergent fixation and allow for 18° of screw angulation

**Square impactation dimple**

Provides rotational control during implantation and functions with a quick release impactor mechanism

**Optimized head to shell ratio**

Allows the use of a 36mm head and liner in a 50mm shell with no compromise to liner thickness

**Self-positioning anti-rotation tabs**


Create an interference fit with the liner to maximize rotational stability while allowing ease of insertion

**BioFoam® cancellous titanium**

Provides immediate fixation and promotes bone apposition

**Minimized shell thickness**

Decreases stiffness, discourages stress-shielding, and allows for optimized liner thickness



# A-Class<sup>®</sup> Highly Cross-Linked Polyethylene

- No oxidation<sup>5</sup>
- Undetectable free radicals<sup>6</sup>
- 92% or greater reduction in wear<sup>7</sup>

**92% or greater reduction in  
WEAR<sup>7</sup>**




■ Conventional Poly

■ A-Class<sup>®</sup> Poly

## Optimized head to shell ratio

The importance of head diameter and its effect on dislocation has been widely recognized in the industry. Clinical data demonstrates that the use of larger femoral heads decrease the risk of dislocation and impingement while increasing range of motion.<sup>8-9</sup>

The Prime Acetabular Cup System offers optimized head to shell ratios and increased jump distance without any compromise to liner thickness. This promotes greater joint stability.



# Robust locking mechanism

The robust locking mechanism coupled with low-profile, self-positioning anti-rotation tabs maximize push-out strength while allowing for easy liner impaction with less than 270 pounds of impaction force required.

- Twelve anti-rotation tabs on the polyethylene liner maximize torque resistance and rotational stability and allow for self-correction during impaction
- Accepts multiple liner configuration options



# Push-out strength

- MicroPort Orthopedics' Prime
- Stryker Trident®
- Smith & Nephew R3™
- Depuy Synthes Pinnacle®

FORCE (IBF)

0

50

100

150

200

250

300

350

400



Pinnacle is a registered trademark of DePuy Synthes, R3 is a trademark of Smith & Nephew, Trident is a registered trademark of Stryker.



# Optimized screw hole locations


Two screw holes are positioned close to the dome hole and one screw hole is positioned close to the periphery, which:

- Allows for divergent screw fixation
- Eliminates the need to adjust inclination or version to avoid vascular injury



# Enhanced Screw Hole Design

- Conical screw hole entry coupled with a spherical screw head allows for poly-axial positioning of the screw anywhere within an 18° cone (9° off-axis)
- Screw holes are recessed to ensure that screws do not sit proud and prevent liner seating at any degree of screw angulation



# Simple, versatile instrumentation


The Prime Acetabular System was designed for compatibility with a variety of surgical approaches, including the following modern, soft tissue-sparing MicroPort surgical approaches:

- Anterior Path®
- SuperPath®
- Path®
- SuperCap®
- Direct Anterior
- Anterolateral
- Posterior



# Modular impaction instrumentation


- Modular impactor handle is compatible with both straight and curved shell and liner impactors
- Shell impactors mate with both a threaded adapter and a quick release adapter that facilitates easy disengagement of the impactor from the shell



# Screw Instrumentation

Screw instrumentation has multiple options to facilitate screw insertion for various surgical approaches and patient anatomies.

- Drill shafts and screw drivers come in standard and long lengths to accommodate standard and high BMI patients, respectively
- Flexible drill shafts and angled screw-holding forceps allow for access through surgical interval in soft tissue-sparing surgical approaches
- Special drills and screw drivers available for use through a portal



# Shell, Liner, and Head Options

The Prime Acetabular Cup System offers various cup and liner options to allow intraoperative flexibility and address a wide variety of patient anatomies.



Available in 46 – 68mm diameters in 2mm increments





STANDARD



LIPPED




+4MM LATERALIZED, 10° FACE-CHANGING

Available in 22, 28, 32, 36, 40, and 44mm diameters



COBALT CHROME METAL



BIOLOX® DELTA® OPTION



BIOLOX® DELTA® CERAMIC



BCH® CERAMIC HEAD



A-CLASS® COBALT CHROME METAL

NOTE: All liners are available in both A-Class® Highly Cross-Linked Polyethylene and E-Class™ Vitamin E Blended Highly Cross-Linked Polyethylene.

**References:**

1. Woodell-May, J.E., Kumar, M. In vitro comparison of cell proliferation on Ti6Al4V and Tantalum Metal. Poster No. 1578 presented at 53rd Annual Meeting of the Orthopaedic Research Society
2. Benchtop data on file at MicroPort Orthopedics.
3. J.D. Bobyn. JBJS [Br] 1999;81-B:907-14.
4. Dyson, E.D., Jackson, C.K. and Whitehouse, W. J., Scanning electron microscope studies of human trabecular bone, Nature, 225, 957-959, 1970.
5. Benchtop data on file at MicroPort Orthopedics.
6. Benchtop data on file at MicroPort Orthopedics.
7. Compared to conventional poly. Benchtop data on file at MicroPort Orthopedics.
8. Benchtop data on file at MicroPort Orthopedics.
9. Bartz, RL et al. The effect of femoral component head size on posterior dislocation of the artificial hip joint. J Bone Joint Surg Am. 2000 Sep; 82 (9): 1300-7.



MicroPort Orthopedics Inc.  
5677 Airline Road  
Arlington, TN USA 38002  
866 872 0211

[microportortho.com](http://microportortho.com)

The CE-Marking of Conformity is applied  
per catalog number and appears on the  
outer package label, if applicable.

Trademarks and Registered marks of  
MicroPort Orthopedics Inc.  
© 2022 MicroPort Orthopedics Inc.  
All Rights Reserved. 012589C JUNE2022