

# *Profemur<sup>®</sup> Renaissance<sup>®</sup>*

*Hip System: Classic and Modular Stems*

Surgical Technique



# Table of Contents

<b>CHAPTER 1</b>	<b>2</b>	<b>DESIGN RATIONALE</b>
<b>CHAPTER 2</b>	<b>3</b>	<b>PRODUCT INFORMATION</b>
	3	Design Features of the Profemur® Renaissance® Classic Stem
	4	Design Features of the Profemur® Renaissance® Modular Stem
	5	Design Features of the Profemur® Modular Neck
	6	General Specifications of the Profemur® Renaissance® Stem
<b>CHAPTER 3</b>	<b>7</b>	<b>PREOPERATIVE PLANNING</b>
<b>CHAPTER 4</b>	<b>8</b>	<b>SURGICAL TECHNIQUE</b>
	8	Femoral Neck Osteotomy
	8	Femoral Canal Preparation
	8	Femoral Reaming
	9	Femoral Broaching
	10	Trial Reduction
	12	Stem Insertion
	13	Implant Assembly
<b>CHAPTER 5</b>	<b>14</b>	<b>TECHNIQUE OVERVIEW</b>
<b>CHAPTER 6</b>	<b>15</b>	<b>IMPLANT REMOVAL</b>
	15	Profemur® Classic Stem Removal
	16	Profemur® Modular Neck and Stem Removal
<b>CHAPTER 7</b>	<b>17</b>	<b>ORDERING INFORMATION</b>
<b>CHAPTER 8</b>	<b>22</b>	<b>INDICATIONS AND WARNINGS</b>

MicroPort Orthopedics recognizes that proper surgical procedures and techniques are the responsibility of the medical professional. The following guidelines are furnished for information purposes only. Each surgeon must evaluate the appropriateness of the procedures based on his or her personal medical training, experience and patient condition. Prior to use of the system, the surgeon should refer to the product package insert for additional warnings, precautions, indications, contraindications and adverse effects. Instructions For Use package inserts are also available by contacting the manufacturer. Contact information can be found on the back of this Surgical Technique and the Instructions For Use package inserts are available on the website listed.

# Chapter 1

## Design Rationale



The Profemur® Renaissance® Hip System is a “ream and broach” stem design that provides immediate stability through fixation of the diaphyseal canal. Reamers prepare the canal diameter to size, and press-fit distal splines provide rotational stability once the stem is implanted. The stem has two medial flare options designed to optimize fit, stability, and load transfer based on the patient anatomy: “standard flare” for femurs with moderate cortical wall thickness, and “reduced flare” for femurs with thick cortical walls (“champagne flute” canals). Secondary fixation is achieved via titanium plasma spray coating to promote bone ingrowth.

The Profemur® Renaissance Hip system has 16 stem sizes to optimize fit for each patient anatomy. Replicating leg length and neck offset is achieved in the Profemur® Renaissance® Classic stem with a combination of short/long and straight/varus necks. For unique patient anatomies requiring even more reconstruction options, the Profemur® Renaissance® Modular stem has 11 neck options of varying lengths and versions. Not only does this modular system optimize stability, range of motion, and leg length, but it also allows for intraoperative versatility to address soft tissue tension as needed<sup>1</sup>.

<sup>1</sup> Traina F, De Clerico M, Biondi F, Pilla F, Tassinari E, Toni A. Sex differences in hip morphology: is stem modularity effective for total hip replacement?; J Bone Joint Surg Am 2009 Nov; 91(6):121-8.

# Chapter 2

## Product Information

### Profemur® Renaissance® Classic Stem Design Features

#### Abbreviated Technique: Ream and Broach

Ream to templated size or cortical chatter

Sequentially broach with reduced flare broaches to corresponding reamer size

If it is not secure, switch to standard flare broach or ream and broach to larger stem size

Implant size and flare corresponding to broach size and flare

#### Ordering Information

Templates	PRNCXR15 (Short Neck Classic)
	PRNLXR15 (Long Neck Classic)
Surgical Technique	008118
Instruments	PRGIKIT1 (Includes trial necks)
	PRZNKIT1 (Broaches and broach handles)
	PRZNKIT2 (Reamers)
	PRCLIMPT (Classic stem impactor)
Implants	PRRCKITS (Short neck stems)
	PRRCKITL (Long neck stems)

#### Sizes

Neck	Short and Long
Reduced Flare	10 - 16
Standard Flare	10 - 18


For additional risk information, please consult the Instructions for Use package insert.

#### Neck Options

Short and long neck lengths, Straight (135° CCD) and Varus 8° (127° CCD) neck angles allowing for multiple head center positions to meet range of anatomical needs


#### Flare

Standard and Reduced flare options designed to maximize metaphyseal fit



#### Proximal Taper

3° proximal anterior/posterior taper designed for proximal fill



#### Driving Platform

Dimple and oval slot designed for unidirectional loading and rotational control during stem insertion, respectively

#### Plasma Spray

Coating thickness provides 0.5mm (0.25mm per side) additional press-fit

#### Surface Roughness

Titanium stem surface has glass-beaded texture

#### Distal Splines

Designed to provide additional 1 mm press-fit (0.5mm per side) for rotational stability



#### Distal Slot

Minimizes stem stiffness to discourage fracture during stem insertion

#### Distal Bullet Tip

Round distal tip designed to reduce the risk of fracture during insertion and minimize point contact after implantation

## Profemur® Renaissance® Modular Stem Design Features

### Abbreviated Technique: Ream and Broach

Ream to templated size or cortical chatter

Sequentially broach with reduced flare broaches to corresponding reamer size

If it is not secure, switch to standard flare broach or ream and broach to larger stem size

Implant size and flare corresponding to broach size and flare


### Ordering Information

Templates	PLS0XR15
Surgical Technique	008118
Instruments	PRGIKIT1 (Includes trial necks) PRZNKIT1 (Broaches and broach handles) PRZNKIT2 (Reamers)
Implants	PLIMKITA (Stems) COCRKITB (Modular necks)

### Sizes


Reduced Flare	10 - 16
Standard Flare	10 - 18

For additional risk information, please consult the Instructions for Use package insert.



#### Driving Platform

Threaded slot designed for rotational control during stem insertion



#### Proximal Taper

3° proximal anterior/posterior taper designed for proximal fill

#### Plasma Spray

Coating thickness provides 0.5mm (0.25mm per side) additional press-fit

#### Flare


Standard and Reduced flare options designed to maximize metaphyseal fit

#### Surface Roughness

Titanium stem surface has glass-beaded texture

#### Distal Splines

Designed to provide additional 1mm press-fit (0.5mm per side) for rotational stability



#### Distal Slot

Minimizes stem stiffness to discourage fracture during stem insertion

#### Distal Bullet Tip

Round distal tip designed to reduce the risk of fracture during insertion and minimize point contact after implantation

## Profemur® Modular Necks Design Features

- Cobalt Chrome Material
- Varus Neck Angle: 127°
- Neutral Neck Angle: 135°
- Valgus Neck Angle: 143°

Dimensional Chart (mm)

Type	Length Options
Neutral	Short Long
Varus/Valgus	Short
Ante/Retro 8 Degree	Short Long
Ante/Retro 15 Degree	Short Long
Ante/Retro - Varus/Valgus 1	Short Long
Ante/Retro - Varus/Valgus 2	Short Long





## Profemur® Renaissance® Stems - General Specifications

### General Specifications

- Stems are made of titanium alloy
- Commercially-pure titanium plasma spray over proximal region (0.5mm/side)
- Medial Stem Length: 125 - 170mm
- M/L Width: 27 – 40mm
- Distal splines are 1mm larger than matching reamer
- A/P Thickness: 14 – 17mm

### Profemur® Renaissance® Hip Stems Dimensional Chart

(Measurements in millimeters)


	Size	Short Neck Classic and Modular		Long Neck Classic and Modular		Stem Measurements			
		Leg Length	Neck Offset	Leg Length	Neck Offset	Med. Length	M/L Width	A/P Thick.	Lat. Length
Reduced	Straight (Neck Angle = 135°)								
	10	31	33	38	41	125	28	14	145
	11	31	35	38	42	129	28	15	148
	12	31	36	39	44	134	30	16	154
	13	32	37	40	45	140	31	17	158
	14	32	37	40	45	145	32	18	163
	15	33	39	41	47	150	33	19	170
Standard	16	33	39	41	47	154	34	21	175
	10	31	33	38	41	125	30	14	145
	11	31	35	38	42	129	30	15	148
	12	31	36	39	44	134	32	16	154
	13	32	37	40	45	140	33	17	158
	14	32	37	40	45	145	34	18	163
	15	33	39	41	47	150	35	19	170
Reduced	16	33	39	41	47	154	36	21	175
	17	33	41	40	49	160	37	22	180
	18	33	42	40	49	170	38	23	190
	Varus 8° (Neck Angle = 127°)								
	10	29	36	35	45	125	28	14	145
	11	27	38	35	46	129	28	15	148
	12	28	40	36	48	134	30	16	154
Standard	13	28	41	37	49	140	31	17	158
	14	28	41	37	49	145	32	18	163
	15	30	43	38	51	150	33	19	170
	16	30	43	38	51	154	34	21	175
	10	29	36	35	45	125	30	14	145
	11	27	38	35	46	129	30	15	148
	12	28	40	36	48	134	32	16	154
Standard	13	28	41	37	49	140	33	17	158
	14	28	41	37	49	145	34	18	163
	15	30	43	38	51	150	35	19	170
	16	30	43	38	51	155	36	21	175
	17	29	45	37	53	160	37	22	180
	18	30	46	37	53	170	38	23	190

Offset and Leg Length are based on +0 head.

### Head Center Adjustment Chart


(Measurements in millimeters)

Head Size	Neck Length Adjustment	OFFSET / LEG LENGTH ADJUSTMENT	
		Straight	Varus 8°
Short	-3.5	-2.5 / -2.5	-2.8 / -2.1
Medium	+0	+0.0 / +0.0	+0.0 / +0.0
Long	+3.5	+2.5 / +2.5	+2.8 / +2.1
X Long	+7	+4.9 / +4.9	+5.6 / +4.2
XX Long	+10.5	+7.4 / +7.4	+8.4 / +6.3



# Chapter 3

## Preoperative Planning



**CAUTION:** Preoperative templating is intended for estimation purposes only. Final component size and position should be determined intraoperatively.

Accurate preoperative templating requires good quality standardized radiographs of the pelvis and operative hip. To determine limb length discrepancy, draw a line across the bottom of the ischium on the A/P view. The distance from this horizontal reference line to each lesser trochanter should then be measured. The difference between each measured side is the leg length discrepancy. If there is any asymmetry of the pelvis or if landmarks are not clear, other means to determine discrepancy should be used.

Determine the femoral head center. Once the center of rotation for the acetabular component has been established, the center of rotation for the femoral head should be determined. Superimpose the femoral stem templates sequentially on the A/P x-ray with the templates positioned neutrally along the longitudinal axis of the femur. Estimate the metaphyseal and diaphyseal fit and anticipated level of implant insertion using the templates. The approximate femoral size and length of the femoral neck cut can be estimated from the templates.

The neck angle and head length which most closely correspond to the patient's femoral head center can be estimated as well. The ideal head will align atop the previously determined center of rotation for the femoral head. In patients with significant deformity of the femoral head, templating can be performed on the opposite hip if necessary.

For soft bone, the implant may seat further than the template indicates. An implant larger than the templated size may be required. For strong, healthy bone, an implant smaller than the templated size may be required.

Each circle represents the center of rotation for a modular short neck with the corresponding head option (Short to XX Long). Each square represents the center of rotation for a modular long neck with the corresponding head option (Short to XX Long). The circles/squares on the AP template of the stem illustrate the impact of choosing an 8° varus/valgus neck relative to the neutral neck position.

*NOTE: AR/VV necks can also affect neck position by 6° varus/valgus.*


The lateral x-ray illustrates the front to back fill of the implant and the position of the implant relative to the femoral anterior bow. If the anterior bow is high, the implant size may be reduced to minimize the risk of fracture. The lateral templates use circles/squares to compare the impact of choosing a neutral neck and necks with 8° or 15° anteversion/retroversion.

Both the A/P and lateral views are needed to illustrate the impact of choosing an AR/VV neck because the combination necks provide multi-dimensional positioning. Each AR/VV neck provides 4° anteversion/retroversion and 6° varus/valgus. The impact of each AR/VV option (1 or 2) depends upon which hip is being considered. Therefore, caution should be used to ensure that the appropriate combination is planned.



# Chapter 4

## Surgical Technique



### Femoral Neck Osteotomy

Using the greater trochanter or lesser trochanter as a reference, resect the neck at a 45° angle to the longitudinal axis of the femur. The Profemur® Neck Resection Guide (P/N PTRG0410) is available to help establish the angle of resection.

### Open the Femoral Canal

Using the Profemur® Box Chisel (P/N PRFS0450), open the femoral canal. The box chisel should be lateralized to ensure a neutral orientation of the implant.

### Starter Reamer

Enter the femoral canal with the Profemur® Initial Reamer (P/N APA04750). Proper reaming depth, measured from the medial resection, corresponds to stem length. Using the Profemur® Renaissance® Proximal Reamer (P/N PSSR1118), clear out the lateral greater trochanter to allow straight access to the femoral canal.

It is important to stay lateral with the box chisel, starter reamer, and proximal reamer. Care should be taken to ensure that the initial reaming track into the femur is in neutral alignment with the femoral axis.


### Femoral Reaming

Starting with the smallest reamer, Profemur® Femoral Reamer Size 8 (P/N PRR00080), sequentially ream the femoral canal to the appropriate depth and diameter based on the templated size and the amount of resistance that the reamer encounters.

Stem length is measured from the upper most plasma line on the medial side (resection line) to the distal tip of the implant. Proper reaming depth, measured off of the medial resection, corresponds to stem length.

### Profemur® Renaissance® Reaming Chart

Stem Size	Reaming Depth
10	125mm
11	130mm
12	135mm
13	140mm
14	145mm
15	150mm
16	155mm
17-18	160mm



Profemur® Neck Resection Guide  
(P/N PTRG0410)



Profemur® Box Chisel  
P/N PRFS0450



Profemur® Initial Reamer  
P/N APA04750



Profemur® Renaissance®  
Proximal Reamer  
P/N PSSR1118



Profemur® Renaissance®  
Reamer Size 8  
P/N PRR00080



### Femoral Broaching

Prepare the femoral canal with the Profemur® Renaissance® Starter Broach (P/N PLRB0009). Staying centered between the anterior and posterior cortices, impact the starter broach until the top of the teeth rests just at or below the level of the neck resection. Insert the broach until it rests 1-2mm below the level of the neck resection.

Attach the broach handle of choice (P/N SLBROHAN is shown) to the smallest reduced flare broach. Begin broaching using a mallet with short, controlled strokes. Sequentially broach to the appropriate size, using the reduced flare Profemur® Renaissance® Broaches (P/N PLRB0010-18).


The correct broach depth is achieved when the base of the polished oval collar rests along the resection. Recognize that the polished collar increases in height as stem size increases. Throughout broaching, continue to apply lateral pressure to ensure neutral alignment of the implant.

Continue broaching until an optimal fit is found. Once the appropriately sized broach is fully seated, leave the broach fully seated in the canal and detach the broach handle to allow for trial reduction.

**CAUTION:** The Profemur® Renaissance® hip system includes both standard and reduced flare broaches and implants. Broaching is performed sequentially, using the reduced flare broaches, to the equivalent reamer size used during the femoral reaming process. If the fit is not secure, sequentially broach to the equivalent reamer size using the standard flare Profemur® Renaissance® Broaches (P/Ns PLSB0010-18).

**NOTE:** The femoral broaches are 0.5mm smaller than the plasma spray of the implant, providing 0.5mm (0.25mm per side) of press-fit proximally.






#### Potential Differences Between Broached and Templated Sizes:

1. The quality of bone plays an integral role in sizing. For soft bone, the broach may seat further than the template indicates. An implant larger than the templated size may be required. Patients with strong, healthy bone might require an implant smaller than the templated size.
2. If a broach smaller than the size templated becomes tight, hard bone at the lateral femoral neck may be pushing the broach into varus. Use the lateral edge of the broach to restore a neutral position. Additional broaching may be necessary.
3. If a broach is going in straight and still becomes tight with sizes smaller than templated, a repetitive in/out broach motion may clear excess medial and lateral bone. If still tight, the stem should be appropriately downsized until metaphyseal bone is engaged.

#### Trial Reduction

Select the appropriate plastic or metal Profemur® Trial Neck (P/Ns APA11102 - APA11154 or APA12102 - APA12154) and either the Profemur® Trial Head (P/Ns APA02121 - APA02154) or the Trial Head/Neck Sleeve combination (P/Ns 41103600 - 41104800 and P/Ns APA0TSS3, APA0TSM0, APA0TSL3 or APA0TSX7) and perform a trial reduction. Metal Profemur® Trial Necks (P/Ns APA12102 - APA12154) are equivalent dimensionally to the plastic Profemur® Trial Necks. Once a well balanced hip has been created with a trial head and trial neck, remove the broach.

For Classic stems, use only the short and long, straight and varus plastic and metal Profemur® Trial Necks (P/Ns APA11102, APA11104, APA11152, APA11154 or APA12102, APA12104, APA12152, APA12154) as they represent the Profemur® Renaissance® Classic stem offerings.



#### Summary of Profemur® Renaissance® Modular Neck Angle Options


The choice of neck anteversion is based on intraoperative assessment of stability. The head/neck combination that allows maximal flexion/internal rotation and extension/external rotation without dislocation should be chosen.

- »» Straight necks create a neutral neck axis.
- »» Varus necks decrease the inclination angle to 127° (neutral position is 135°); the femoral head shifts medially and inferiorly; leg length is shortened; offset is increased.
- »» Valgus necks increase the inclination angle to 143°; the femoral head shifts laterally and superiorly; leg length is increased; offset is decreased.
- »» Anteverted necks shift the femoral head anteriorly relative to the stem by 8° or 15°.

- »» Retroverted necks shift the femoral head posteriorly relative to the stem by 8° or 15°. Retroverted necks prove useful in hips with excess femoral anteversion such as DDH.
- »» AR/VV necks combine anteversion/retroversion and varus/valgus necks to offer a broad range of multi-dimensional head positions. Each AR/VV neck provides 4° of A/R and 6° of V/V.

#### Summary of Profemur® Renaissance® Classic Neck Angle Options

- »» Straight (135°) necks create a neutral neck axis.
- »» Varus 8° necks decrease the inclination angle to 127°; the femoral head shifts medially and inferiorly; leg length is shortened; offset is increased.



## Stem Insertion

### Modular Stem Option 1

Connect the Profemur® Threaded Inline Stem Inserter (P/N PRFS1461) into the small threaded hole on the proximal face of the femoral implant, and insert the stem into the femoral canal with hand pressure to align proper version. Use a mallet with short, controlled strokes to seat the implant. Typically, the final implant is seated with the base of the polished collar aligned with the neck resection.

### Modular Stem Option 2

Insert the femoral implant into the femoral canal and align the version before placing the Profemur® Screwdriver Inserter (P/N PRFS0460) into the impaction slot on the proximal face. Use a mallet with short, controlled strokes to seat the implant. Typically, the final implant is seated with the base of the polished collar aligned with the neck resection.

### Classic Stem Option 1

The Profemur® Classic Stem Inserter (P/N PRCLIMPT) is available to provide rotational control during impaction.

Place the tip of this impactor into the oval slot on the proximal face and, with a mallet, fully seat the implant using short, controlled strokes. Typically, the final implant is seated with the base of the polished collar aligned with the neck resection.

### Classic Stem Option 2

Place the tip of the Final Stem Impactor (P/N PPF60200) into the dimple on the proximal face and, with a mallet, fully seat the implant using short, controlled strokes. Typically, the final implant is seated with the base of the polished collar aligned with the neck resection.

### Final Trial Reduction

Perform a final trial reduction using the plastic necks and trial heads to reconfirm stability, range of motion and leg length.

**CAUTION:** Do not use metal trial necks with the implant. Metal trial necks are only to be used with broaches since they may damage the taper. Only plastic trial necks (available in PRGKIT1) should be used for trial reductions with the implant.




### Implant Assembly

To properly assemble and impact a Profemur® modular neck, the following procedure is recommended:

**STEP A.** Suction any fluid from the stem implant pocket. Ensure that both the stem and neck are clean and dry prior to assembly.

**STEP B.** Insert the oval end of the appropriate femoral neck implant into the femoral stem pocket.


**STEP C.** Position the leg such that the knee is supported by an assistant on the opposite side of the table. By resting the patient's knee against the mid-section of the assistant, this will provide counter-force against the mallet blows to ensure the impaction load transfer to the neck junction.



**STEP D.** Affix the femoral head to the neck. Using the head impactor instrument, strike the impactor with **three very firm blows** with a mallet to securely fix the head to the neck and stem.

*NOTE: Place a femoral head impactor (such as P/N 4400F10000 or PPR67200) on a ceramic femoral head, and align the impactor with the femoral neck axis of the stem implant. With a moderate top of the hammer in the axial direction, firmly impact the ceramic femoral head on the stem taper until it is fully seated.*

*NOTE: If using a Profemur® Renaissance® Classic stem, affix the femoral head to the stem and impact as instructed.*







# Chapter 5

## Technique Overview


1. X-ray



2. Femoral Neck Osteotomy



3. Femoral Canal Preparation




4. Femoral Reaming





5. Femoral Broaching




6. Trial Reduction



7. Stem Insertion



8. Implant Assembly



# Chapter 6

## Implant Removal

### Profemur® Renaissance® Classic Stem Removal



*If the removal of the implant is required due to revision or failure of the device, the surgeon should contact the manufacturer using the contact information located on the back cover of this surgical technique to receive instructions for returning the explanted device to the manufacturer for investigation.*

#### Classic Stem Removal

Should the removal of a Profemur® Classic stem become necessary, the Universal Stem Extractor (P/N 4700SE05) and the corresponding Slap Hammer (P/N 4700SH0000) can be utilized. Thread the stem extractor onto the threaded end of the slap hammer. With the femoral head removed, position the stem extractor across the flats on the sides of the femoral neck, and remove the stem using repetitive upward blows delivered by the slap hammer.

*NOTE: Classic stem extraction instruments must be ordered separately.*



Universal Stem Extractor  
P/N 4700SE05



Slap Hammer  
P/N 4700SH0000

## Profemur® Renaissance® Neck and Stem Removal

### Modular Neck Removal

Should it be necessary, a Profemur® hip stem can be removed in the following manner.

**STEP A.** The femoral head is removed by placing an osteotome or bone punch on the underside of the femoral head and applying mallet blows upward until the femoral head is removed. With the femoral head removed, thread the 12/14 Adapter (P/N APA00003) over the round taper end of the modular neck.




**STEP B.** Place the Head/Neck Extractor (P/N APA00001) over the 12/14 adaptor and modular neck, and hand tighten the hex end of the shaft until the base of the extractor rests on the stem, while the fork of the extractor rests under the rim of the adaptor.



**CAUTION:** The base of the extractor must rest on the top surface of the stem's modular neck pocket, and not on the resected bone.



**STEP C.** Attach the Spanner Handle (P/N APA00005) to the hex end of the extractor and rotate clockwise until the neck taper disengages. The Tommy Bar (P/N APA00006) can be inserted into the end of the spanner handle for even greater leverage.



Please note that these instruments are designed for the purpose of removing a neck during the primary surgery.

These instruments may or may not be able to provide the force necessary to disengage a connection between components that have been implanted for a longer period of time. In revision cases, removal and replacement of only the modular neck is contraindicated.

### Modular Stem Removal

The thread at the base of a Profemur® stem's modular neck pocket can now be accessed to remove the stem. Insert the Femoral Stem Extractor (P/N PPR67688) into the modular neck pocket and tighten the threaded shaft by hand, firmly seating the shaft via the use of the Hex Screwdriver (P/N PP275400). Using the slide hammer, create extraction forces onto the underside of the shaft using repetitive upward blows to remove the stem. If bone on-growth exists, it may be necessary to use osteotomes in order to first disengage the stem/bone interface.



12/14 Adapter  
P/N APA00003



Head/Neck Extractor  
P/N APA00001



Spanner Handle  
P/N APA00005



Tommy Bar  
P/N APA00006



Hex Screwdriver  
P/N PP275400



Femoral Stem Extractor  
P/N PPR67688

# Chapter 6

## Ordering Information



### Profemur® Renaissance® Classic Stems - Short Necks

#### PRRCKITS

Catalog No.	Flare	Neck Angle	Size
PLSES010	Standard	Varus 8°	10
PLSES011	Standard	Varus 8°	11
PLSES012	Standard	Varus 8°	12
PLSES013	Standard	Varus 8°	13
PLSES014	Standard	Varus 8°	14
PLSES015	Standard	Varus 8°	15
PLSES016	Standard	Varus 8°	16
PLSES017	Standard	Varus 8°	17
PLSES018	Standard	Varus 8°	18
PLSSS010	Standard	Straight	10
PLSSS011	Standard	Straight	11
PLSSS012	Standard	Straight	12
PLSSS013	Standard	Straight	13
PLSSS014	Standard	Straight	14
PLSSS015	Standard	Straight	15
PLSSS016	Standard	Straight	16
PLSSS017	Standard	Straight	17
PLSSS018	Standard	Straight	18
PLSER010	Reduced	Varus 8°	10
PLSER011	Reduced	Varus 8°	11
PLSER012	Reduced	Varus 8°	12
PLSER013	Reduced	Varus 8°	13
PLSER014	Reduced	Varus 8°	14
PLSER015	Reduced	Varus 8°	15
PLSER016	Reduced	Varus 8°	16
PLSSR010	Reduced	Straight	10
PLSSR011	Reduced	Straight	11
PLSSR012	Reduced	Straight	12
PLSSR013	Reduced	Straight	13
PLSSR014	Reduced	Straight	14
PLSSR015	Reduced	Straight	15
PLSSR016	Reduced	Straight	16



### Profemur® Renaissance® Classic Stems - Long Necks

#### PRRCKITL

Catalog No.	Flare	Neck Angle	Size
PLLES010	Standard	Varus 8°	10
PLLES011	Standard	Varus 8°	11
PLLES012	Standard	Varus 8°	12
PLLES013	Standard	Varus 8°	13
PLLES014	Standard	Varus 8°	14
PLLES015	Standard	Varus 8°	15
PLLES016	Standard	Varus 8°	16
PLLES017	Standard	Varus 8°	17
PLLES018	Standard	Varus 8°	18
PLLSS010	Standard	Straight	10
PLLSS011	Standard	Straight	11
PLLSS012	Standard	Straight	12
PLLSS013	Standard	Straight	13
PLLSS014	Standard	Straight	14
PLLSS015	Standard	Straight	15
PLLSS016	Standard	Straight	16
PLLSS017	Standard	Straight	17
PLLSS018	Standard	Straight	18
PLLER010	Reduced	Varus 8°	10
PLLER011	Reduced	Varus 8°	11
PLLER012	Reduced	Varus 8°	12
PLLER013	Reduced	Varus 8°	13
PLLER014	Reduced	Varus 8°	14
PLLER015	Reduced	Varus 8°	15
PLLER016	Reduced	Varus 8°	16
PLLSR010	Reduced	Straight	10
PLLSR011	Reduced	Straight	11
PLLSR012	Reduced	Straight	12
PLLSR013	Reduced	Straight	13
PLLSR014	Reduced	Straight	14
PLLSR015	Reduced	Straight	15
PLLSR016	Reduced	Straight	16



## Profemur® Renaissance® Stems

### PLIMKITA

Catalog No.	Description
PLS0S410	Profemur® Renaissance® Standard Flare Size 10
PLS0S410	Profemur® Renaissance® Standard Flare Size 11
PLS0S410	Profemur® Renaissance® Standard Flare Size 12
PLS0S410	Profemur® Renaissance® Standard Flare Size 13
PLS0S410	Profemur® Renaissance® Standard Flare Size 14
PLS0S410	Profemur® Renaissance® Standard Flare Size 15
PLS0S410	Profemur® Renaissance® Standard Flare Size 16
PLS0S410	Profemur® Renaissance® Standard Flare Size 17
PLS0S410	Profemur® Renaissance® Standard Flare Size 18
PLSSS010	Profemur® Renaissance® Reduced Flare Size 10
PLSSS011	Profemur® Renaissance® Reduced Flare Size 11
PLSSS012	Profemur® Renaissance® Reduced Flare Size 12
PLSSS013	Profemur® Renaissance® Reduced Flare Size 13
PLSSS014	Profemur® Renaissance® Reduced Flare Size 14
PLSSS015	Profemur® Renaissance® Reduced Flare Size 15
PLSSS016	Profemur® Renaissance® Reduced Flare Size 16






## Profemur® Modular Necks

### COCRKITB

Catalog No.	Description
PHAC1202	Straight Short CoCr
PHAC1204	Straight Long CoCr
PHAC1212	Ante/Retro - Varus/Valgus 2 Short CoCr
PHAC1214	Ante/Retro - Varus/Valgus 2 Long CoCr
PHAC1222	Ante/Retro - Varus/Valgus 1 Short CoCr
PHAC1224	Ante/Retro - Varus/Valgus 1 Long CoCr
PHAC1232	Ante/Retro 8° Short CoCr
PHAC1234	Ante/Retro 8° Long CoCr
PHAC1242	Ante/Retro 15° Short CoCr
PHAC1244	Ante/Retro 15° Long CoCr
PHAC1252	Varus/Valgus 8° Short CoCr

# Profemur® Standard Instrument Kit

PRGIKIT1



Catalog No.	Description
20070050	Modular Neck Inserter
4400FI0000	Femoral Head Impactor
APA00001	Head/Neck Extractor
APA00003	Head/Neck Extractor Adaptor 12/14
APA00005	Head/Neck Extractor Cardan Spanner Hex
APA00006	Head/Neck Extractor Tommy Bar
APA02121	Femoral Trial Head 28mm Short (-3.5mm)
APA02122	Femoral Trial Head 28mm Medium (+0mm)
APA02123	Femoral Trial Head 28mm Long (+3.5mm)
APA02124	Femoral Trial Head 28mm XLong (+7mm)
APA02125	Femoral Trial Head 28mm XXLong (+10.5mm)
APA02142	Femoral Trial Head 36mm Short (-3.5mm)
APA02144	Femoral Trial Head 36mm Medium (+0mm)
APA02146	Femoral Trial Head 36mm Long (+3.5mm)
APA02148	Femoral Trial Head 36mm XLong (+7mm)
APA02151	Femoral Trial Head 32mm Short (-3.5mm)
APA02152	Femoral Trial Head 32mm Medium (+0mm)
APA02153	Femoral Trial Head 32mm Long (+3.5mm)
APA02154	Femoral Trial Head 32mm XLong (+7mm)
APA04241	Profemur® MIS Broach Handle (Qty 2)
APA04244	Broach Handle Alignment Guide Rod (Qty 2)
APA04750	Profemur® Starter Reamer
APA11102	Profemur® Short Straight Plastic Trial Neck
APA11104	Profemur® Long Straight Plastic Trial Neck
APA11112	Profemur® Short A/R Var/Val 1 Plastic Trial Neck
APA11114	Profemur® Long A/R Var/Val 1 Plastic Trial Neck
APA11122	Profemur® Short A/R Var/Val 2 Plastic Trial Neck
APA11124	Profemur® Long A/R Var/Val 2 Plastic Trial Neck
APA11132	Profemur® Short A/R 8° Plastic Trial Neck
APA11134	Profemur® Long A/R 8° Plastic Trial Neck
APA11142	Profemur® Short A/R 15° Plastic Trial Neck
APA11144	Profemur® Long A/R 15° Plastic Trial Neck
APA11152	Profemur® Short Var/Val 8° Plastic Trial Neck
APA11154	Profemur® Long Var/Val 8° Plastic Trial Neck
K0001016	Quick Disconnect T-Handle
PP275400	Hex Screwdriver
PPR67688	Slap Hammer Stem Extractor
PRFS0450	Profemur® Box Chisel
PRFS0460	Profemur® Screwdriver Inserter
PRFS0462	Profemur® Broach Extraction Shaft
PRFS0463	Profemur® Tissue Protecting Sleeve
PRFS1461	Profemur® Threaded In-Line Stem Inserter





## Instruments

### PRZNKIT1

Catalog No.	Description
PLRB0009	Profemur® Renaissance® Starter Broach
PLSB0010	Profemur® Renaissance® Broach Standard 10
PLSB0011	Profemur® Renaissance® Broach Standard 11
PLSB0012	Profemur® Renaissance® Broach Standard 12
PLSB0013	Profemur® Renaissance® Broach Standard 13
PLSB0014	Profemur® Renaissance® Broach Standard 14
PLSB0015	Profemur® Renaissance® Broach Standard 15
PLSB0016	Profemur® Renaissance® Broach Standard 16
PLSB0017	Profemur® Renaissance® Broach Standard 17
PLSB0018	Profemur® Renaissance® Broach Standard 18
PLRB0010	Profemur® Renaissance® Broach Reduced 10
PLRB0011	Profemur® Renaissance® Broach Reduced 11
PLRB0012	Profemur® Renaissance® Broach Reduced 12
PLRB0013	Profemur® Renaissance® Broach Reduced 13
PLRB0014	Profemur® Renaissance® Broach Reduced 14
PLRB0015	Profemur® Renaissance® Broach Reduced 15
PLRB0016	Profemur® Renaissance® Broach Reduced 16
PLRB0017	Profemur® Renaissance® Broach Reduced 17
PLRB0018	Profemur® Renaissance® Broach Reduced 18
PTRG0410	Profemur® Renaissance® Neck Resection Guide
PLRB0016	Profemur® Slotted Broach Impactor Handle (Qty 2)

## Instruments

### PRZNKIT2

Catalog No.	Description
PSSR1118	Profemur® Renaissance® Proximal Reamer
PRR00080	Profemur® Femoral Reamer 8
PRR00085	Profemur® Femoral Reamer 8.5
PRR00090	Profemur® Femoral Reamer 9
PRR00095	Profemur® Femoral Reamer 9.5
PRR00100	Profemur® Femoral Reamer 10
PRR00105	Profemur® Femoral Reamer 10.5
PRR00110	Profemur® Femoral Reamer 11
PRR00115	Profemur® Femoral Reamer 11.5
PRR00120	Profemur® Femoral Reamer 12
PRR00125	Profemur® Femoral Reamer 12.5
PRR00130	Profemur® Femoral Reamer 13
PRR00135	Profemur® Femoral Reamer 13.5
PRR00140	Profemur® Femoral Reamer 14
PRR00145	Profemur® Femoral Reamer 14.5
PRR00150	Profemur® Femoral Reamer 15
PRR00155	Profemur® Femoral Reamer 15.5
PRR00160	Profemur® Femoral Reamer 16
PRR00165	Profemur® Femoral Reamer 16.5
PRR00170	Profemur® Femoral Reamer 17
PRR00175	Profemur® Femoral Reamer 17.5
PRR00180	Profemur® Femoral Reamer 18





## Stem Impactors

Catalog No.	Description
PRCLIMPT	Profemur® Classic Stem Inserter
PRFS1461	Profemur® Threaded Inline Stem Inserter
PPF60200	Final Stem Impactor
PRFS0460	Profemur® Screwdriver Inserter




## Woodpecker Broaching System

4251KT10

Catalog No.	Description
APA00930	Woodpecker
APA00931	Hose Attachment
APA00943	Broach Handle

## Profemur® Renaissance® X-Ray Templates

Catalog No.	Description
PLS0XR15	Profemur® Renaissance® Modular X-ray Template 15% Magnification
PRNCXR15	Profemur® Renaissance® Classic Short Neck X-ray Template 15% Magnification
PRNLXR15	Profemur® Renaissance® Classic Long Neck X-ray Template 15% Magnification



## Broach Handle

SPBHKIT1

Catalog No.	Description
INLNBRHN	Inline Broach Handle (Qty 2)

# Chapter 7

## Indications and Warnings

### Intended Use

MicroPort total hip systems are intended for use in total hip arthroplasty for reduction or relief of pain and/or improved hip function in skeletally mature patients.

### Indications for Use

- 1) non-inflammatory degenerative joint disease such as osteoarthritis, avascular necrosis, ankylosis, protrusio acetabuli, and painful hip dysplasia;
- 2) inflammatory degenerative joint disease such as rheumatoid arthritis;
- 3) correction of functional deformity; and,
- 4) revision procedures where other treatments or devices have failed

Rough grit blast surfaces and the titanium plasma spray coatings applied to implant surfaces are intended for uncemented arthroplasty.

### Contraindications

Patients should be warned of these contraindications.

Contraindications include:

- 1) overt infection;
- 2) distant foci of infections (which may cause hematogenous spread to the implant site);
- 3) rapid disease progression as manifested by joint destruction or bone absorption apparent on roentgenogram;
- 4) skeletally immature patients (patient is less than 21 years of age at the time of surgery);
- 5) cases where there is inadequate neuromuscular status (e.g., prior paralysis, fusion and/or inadequate abductor strength), poor bone stock, poor skin coverage around the joint which would make the procedure unjustifiable;
- 6) neuropathic joints;
- 7) hepatitis or HIV infection;
- 8) neurological or musculoskeletal disease that may adversely affect gait or weight-bearing.

### Product-Specific Warnings and Precautions

Do not attempt to seat the implant beyond the envelope of femoral bone preparation. Forcing to seat the implant beyond the prepared femoral bone may increase the chance of bone fracture. In some cases, a portion of the proximal body with or without coating may be visible above the proximal resection level.

The smaller sized femoral implants are intended for patients with narrower intramedullary femoral canals. The geometry of these implants is reduced to accommodate the anatomy of the narrower intramedullary femoral canal, which also decreases the fatigue-strength and load-bearing characteristics of the implant.

### Other Modular Components (Femoral Head and Stems, Modular Necks and Proximal Body)

Scratching of femoral heads, modular necks and proximal and distal stem tapers should be avoided. Repeated assembly and disassembly of these components could compromise the locking action of the taper joint. Prior to assembly, surgical debris must be cleaned from the interior of the female seat of the proximal body

to ensure proper locking. Ensure components are firmly seated to prevent disassociation. The femoral head, neck taper of the femoral component, modular neck tapers, body taper, female seat of the proximal body **must** be clean and dry before assembly. Do not resterilize femoral prostheses with ceramic femoral heads seated on the stem. Please refer to the product package insert for specific warnings and precautions regarding ceramic femoral heads.

Stems and modular necks with the MicroPort 12/14 SLT Taper should only be used in combination with femoral heads with the MicroPort 12/14 SLT Taper. Cobalt chrome femoral heads with the MicroPort 12/14 SLT Taper are designed for use with cobalt-chromium-molybdenum, titanium alloy and ISO 5832-9 stainless steel (not available in the U.S. or Canada) femoral components with the MicroPort 12/14 SLT Taper.

The neck/body component or neck/femoral stem should be changed only when clinically necessary. Refer to proper neck extraction technique in the surgical technique.

### Modular Necks

- Cobalt Chrome Modular Necks are not for use with the following devices:
  - o Alumina (Bilox Forte) "Ceramic Femoral Head" (size 28mm Long)
- In the United States, **ONLY** the following stems are cleared to be used with both options of titanium or cobalt chrome modular necks:
  - o Profemur® R Stem
  - o Profemur® Z Grit Blast Stem
  - o Profemur® Renaissance® Stem
  - o Profemur® LX Revision Stem
  - o Profemur® TL Stem

All other stems are cleared for use with cobalt chrome modular necks only.

The potential long-term biological effects of metal wear debris and metal ion production are not known. Questions regarding carcinogenicity have been raised in literature; no studies have conclusive evidence that metal wear debris or metal ions are carcinogenic.

**NEVER** combine modular or hard bearing components made by different manufacturers.

Ceramic femoral heads should not be placed on scratched or previously assembled metal tapers as this may lead to a ceramic fracture.

### IMPORTANT

Prior to use of the system, the surgeon should refer to the product package insert for additional warnings, precautions, indications, contraindications and adverse effects. Instructions For Use package inserts are also available by contacting the manufacturer. Contact information can be found on the back of this Surgical Technique and the Instructions For Use package inserts are available on the website listed.




*Full Function. Faster®*



MicroPort Orthopedics Inc.  
5677 Airline Road  
Arlington, TN USA 38002  
866 872 0211

[ortho.microport.com](http://ortho.microport.com)



MicroPort Orthopedics BV  
Hoogoorddreef 5  
1101 BA Amsterdam  
The Netherlands  
+31 20 545 01 00

The CE-Marking of Conformity is applied per catalog number and appears on the outer package label, if applicable.

Trademarks and Registered marks of MicroPort Orthopedics Inc.  
© 2017 MicroPort Orthopedics Inc. All Rights Reserved. 008118 DEC2017